

In This Issue

[From the Editor's Basket](#)

[Thermal Conductivity Detectors](#)

[Announcements](#)

[OEM TCD Supplier](#)

[TCD Troubleshooting - Part 1](#)

[Replacement Thermal Conductivity Detectors](#)

["Basic Gas Chromatography: Techniques in Analytical Chemistry" Textbook](#)

[Trade in your Model 69-100 GC](#)



GOW-MAC Instrument Co.
World Headquarters, Bethlehem, PA U.S.A.

From the Editor's Basket

In this issue of the "Hotwire", I'll be discussing the thermal conductivity detector (TCD). The TCD is one of the more widely used detectors for gas chromatography and gas analysis because it is simple in construction, rugged, versatile, sensitive, relatively linear over a wide range, nondestructive to the sample, and inexpensive. It is because of these qualities that it continues to be widely used even though more sensitive and specialized detectors have been developed.

GOW-MAC has been manufacturing TCDs for well over 60 years now. We make more, and have made more, thermal conductivity detectors than anyone else in the world. Chances are we supplied the one that's in your GC or gas analyzer (or at least the filaments) at this very moment, regardless if it's a GOW-MAC instrument or not.

Contact us to get the answers regarding thermal conductivity detectors. Whether you need to replace one in an existing instrument or need to select one for a special application, we'll make sure it's right for you.


We hope you enjoy this latest issue of the "Hotwire". Your comments and suggestions are always welcome.



The Thermal Conductivity Detector (TCD)

The [Thermal Conductivity Detector \(TCD\)](#) can contain two types of transducers 1) hot-wire filaments or 2) thermistors. The vast majority of TCD's in use today are of the hot-wire type. Historically they were used in the thermal conductivity cell in the 1840's when the effect of hydrogen was observed on a glowing wire and noted in the proceedings of the Royal Irish Academy. Thermistors were developed at Bell Telephone Laboratories in the mid 1940's.

The TCD also called katharometer, is a metal block with paths drilled for gas flow. The detector element is located in this path or cavity and is heated by an electric current. The loss of heat from the element to the cell wall is determined by the ability of the gas in the cavity to conduct heat, i.e., the thermal conductivity of the gaseous medium in the cavity. The elements are connected in a Wheatstone bridge circuit, one exposed to a reference gas, the other to a sample gas. Other designs have developed which include four elements, eight elements, novel flow patterns, heat sinks and completion resistors.



continued on page 2

OEM TCD Supplier

GOW-MAC makes more thermal conductivity detectors than anyone in the world. Our Engineering Department will work on a proprietary basis with you to design a detector for your specific needs. Our detectors have gone into space as well as have been incorporated into CAT Scan machines.



A TCD is a rugged, small, fast response sensor that is inexpensive to purchase and operate. The TCD is concentration dependent and universal. If you think your application or instrument requires a TCD, call us at (610) 954-9000 and discuss your idea in confidence.


TCD TROUBLE SHOOTING GUIDE - Part 1

| TROUBLE | PROBABLE CAUSE | CHECKS AND/OR REMEDY |
|-------------------------------|--|--|
| 1. No signal | a) Detector or power switch off. | Make sure carrier is flowing and turned on. Set bridge current to desired setting. |
| | b) Recording device improperly connected. | Check connections. |
| | c) Detector filaments burned out. | Replace elements. |
| | d) Open circuit on detector cell. | Check for broken lead wires, loose terminal screws and broken filament(s). |
| 2. Low cell current capacity. | a) Power Supply voltage inadequate (likely to occur with change from lower resistance element to higher resistance element). | Increase power supply. |
| 3. Recorder can't be zeroed. | a) Excessive bridge current. | Reduce bridge current. |
| | b) Detector contaminated. | Clean cell cavities and elements. Replace elements if required. |

continued on page 4

Thermistors (cont.)

Under certain conditions thermistor elements can provide higher sensitivity than metal filaments, but the conditions are so restrictive that thermistors cannot effectively perform as wide-range, general purpose detectors. Hotwire elements provide higher overall sensitivity, superior stability, longer operating life, fewer anomalies of operation and lower cost replacement.



GOW-MAC Makes Replacement Thermal Conductivity Detectors

GOW-MAC has been supplying TCD's and filaments to Gas Chromatography manufacturers for over 40 years. Manufacturer's cells we repair or replace include:

Agilent (Hewlett Packard, F & M)
 Beckman
 Bendix
 Burrell
 Carlo Erba
 Chrompak (Packard)
 Finnigan (Tremetrics, Tracor)
 Fisher Scientific (Victoreen)
 Intersmat
 Nuclear Chicago (Barber Coleman)
 Perkin Elmer
 Philips Electronic
 Siemens
 Toyo
 Varian

If the manufacturer of your TCD is not listed, contact us with the manufacturer and model number and we'll determine if [replacement filaments](#) or a replacement TCD is available.

U.S.A. / Int'l Headquarters & Plant

GOW-MAC Instrument Co.
 277 Brodhead Road
 Bethlehem, PA 18017-8600 U.S.A.
 Tel: (610) 954-9000
 Fax: (610) 954-0599
 E-mail: repairs@gow-mac.com

European Office & Plant

GOW-MAC Instrument Co., Ltd.
 Bay K, 14a Industrial Estate
 Shannon, Cty. Clare, Ireland
 Tel: +353-61-471632
 Fax: +353-61-471042
 E-mail: sales@gow-mac.ie

| <u>TROUBLE</u> | <u>PROBABLE CAUSE</u> | <u>CHECKS AND/OR REMEDY</u> |
|-------------------------------------|---|--|
| 3. Recorder can't be zeroed (cont.) | c) Loose or corroded electrical connections. | Check connections. |
| | d) Detector elements oxidized bridge out of balance. | Replace all elements. |
| 4. Drift. | a) Change in flow sample and/or reference gas pressure. | Check flow gauges. Creeping regulator or reducing valve. |
| | b) Warm-up period too short. | Allow sufficient time for purge and temperature equilibration. |
| | c) Changes in ambient temperature. | Protect instrument from drafts, direct sunlight, or nearby sources of hot or cool air. |
| | d) Carrier gas flow leaks . | Tighten all fittings so they are leak free. |
| | e) Detector contaminated. | Clean detector cell. |
| | f) Contamination in column; Column bleed. | Recondition column. |
| | g) Filament aging. | Replace filaments. |
| | h) Decompression chill by the reference gas. | Install heater or buffer volume. |

Find you need to quickly get up to speed on the principals and practice of modern gas chromatography? Then you'll need this:

"Basic Gas Chromatography: Techniques in Analytical Chemistry" by Drs. McNair and Miller, 1998.

The book is designed to serve as a primer/working reference for bench chemists and as a textbook for upper-level undergraduate and graduate students. Written by two well-known practitioners and educators in GC, this book offers thorough coverage of the basics of GC.

Click [here](#) to read more.


Visa®, MasterCard®, and American Express® accepted.

Attention Chemistry Departments!

Trade in your old GOW-MAC Model 69-100 GC

When you purchased your GOW-MAC Model 69-100 Gas Chromatograph in 1969, it cost you \$595. You can now trade the 69-100 in and receive a credit of **\$300.00** towards the purchase of a new [Series 350 or Series 400](#) Thermal Conductivity Detector Gas Chromatograph and [Chrom Perfect® Data Acquisition Software Package](#).

Talk about value – where else can you purchase a GC, use it for 30 years and get a credit of 50% of the original sales price back. The GC has cost you less than \$10.00 per year. Your new Series 350 or 400 will give you another 30+ years of service.



Model 69-100 Gas Chromatograph

So, search those cupboards, closets and deep dark recesses. If you find a GOW-MAC 69-100 GC give us a call at **(610) 954-9000** and ask for Ken.